

Hy Con Check

Kits for Environmental Surface Sampling

Hy Con Check line is a range of innovative solutions including environmental sampling sponges, wipes, and ready-to-use sampling kits, designed for environmental surface sampling of work areas, equipment, animal carcasses, and any other place where testing for Listeria, Salmonella, E. coli and other food-borne pathogens is required. The special materials are free of bactericides and has been tested to be non-inhibitory.

Sterile Sampling Materials

Meccontti provides a wide range of irradiated, non-inhibitory materials which allows to choose the most suited tools to collect samples from the surface...

Pre-moistened or dry formats available to fit many applications













For many different purposes

Animal Carcasses Sampling

Environmental surface sampling

Food-borne pathogenes testing



Code	Product	Product description	Package	
SK001.0010	Hy Con Check W Dry sterile sampling wipes	Content of 1 set- irradiated viscose cloth (25x40cm)to collect environmental samples, glove and the sample bag with a pocket for the description of samples. Sampling wipes can be used with our neutralizing buffer, peptone water, or any other medium used by the client. Because of size, cloths are useful for taking samples from the larger surfaces, production equipment or parts. With peptone water may be used in poultry houses or other farm buildings, for collecting samples against Salmonella from drinking fountains or other elements.	1 bag -10 sets	
SK003.0005	Hy Con Check S Dry sterile sampling sponges	Content of 1 set- irradiated cellulose sponge (size 82x61mm) to collect environmental samples, glove and sample bag with a pocket for the description of samples.	1 bag - 5 sets	
SK004.0005	Hy Con Check S/BPW Sterile sampling sponges pre-moistened with BPW	Content of 1 set- irradiated, cellulose sponge (size 82x61mm) to collect environmental samples (pre-moistened with 10ml of Buffered Peptone Water), glove and a sample bag with a pocket for the description of samples . Sampling sponges are pre-moistened with Buffered Peptone Water. Recommended for Red Meat Carcass Sampling.	1 bag -5 sets	
SK005.0010	Hy Con Check S/BPW Kit Kit of dry sterile sampling sponges and BPW for sponges moistening	Content of 1 set- irradiated cellulose sponge (size 82x61mm) to collect environmental samples, glove, sample bag with a pocket for the description of samples and 1 vial of 10ml Buffered Peptone Water .	1 kit 1 bag of 10 sets SK003 + 10 vials of 10ml BPW	
SK006.0010	Hy Con Check W/BPW Kit Kit of dry sterile sampling wipes and BPW for wipes moistening	Content of 1 set- irradiated viscose cloth (25x40cm) to collect environmental samples, glove, sample bag with a pocket for the description of samples and 1 vial of 10ml Buffered Peptone Water.	1 kit 1 bag of 10 sets SK001 + 10 vials of 10ml BPW	
SK008.0010	Hy Con Check W/NPB Kit Kit of dry sterile sampling wipes and NPB for wipes moistening	Content of 1 set- irradiated viscose cloth (25x40cm) to collect environmental samples, glove, sample bag with a pocket for the description of samples and 1 vial of 10ml Neutralizing Peptone Buffer+LTTh.	1 kit 1 bag of 10 sets SK001 + 10 vials of 10ml NPB)	
SK009.0010	Hy Con Check S/NPB Kit Kit of dry sterile sampling sponges and NPB for sponges moistening	Content of 1 set- irradiated cellulose sponge (size 82x61mm) to collect environmental samples, glove, sample bag with a pocket for the description of samples and 1 vial of 10ml Neutralizing Peptone Buffer+LTTh.	1 kit 1 bag of 10 sets SK003 + 10 vials of 10ml NPB	
SK010.0005	Hy Con Check S Dry sterile sampling sponges	Content of 1 set- irradiated cellulose sponge (size 82x61mm) to collect environmental samples, zip bag with a sticker for the sample description	1 bag -5 sets	
SK012.0005	Hy Con Check S/BNaCl Sterile sampling sponges pre-moistened with BNaCl	Content of 1 set- irradiated cellulose sponge (size 82x61mm) to collect environmental samples (pre-moistened with 10ml of (NaCL - Peptone solution, buffered), glove and a sample bag with a pocket (bag size: 180x250mm) for the description of samples.	1 bag - 5 sets	
SK015.0005	Hy Con Check S/NPB Sterile sampling sponges pre-moistened with NPB	Content of 1 set- irradiated cellulose sponge (size 82x61mm) to collect environmental samples (pre-moistened with 10ml of NPB - Neutralizing Peptone Buffer) glove and a sample bag with a pocket (bag size:185X250mm) for the description of samples. Recommended for environmental samples where disinfectants were used. Neutralizing Buffer will neutralize most of them and provide the end user with total testing assurance.	1 bag - 5 sets	
SK030.0001	Hy Con Check template	Sterile, flexible, individually packed cardboard templates for easy sampling. Recommended for use with Mecconti Hy Con Check sampling kits. Available in size:10 x 10cm. 	1 template min order 10 pcs	
MCBPW10	Sampling media Buffered Peptone Water	In many cases a surface will contain traces of disinfectants and/or sanitizers, but it also possible that it will contain a small number of bacteria too. In such cases, to avoid false-negative results we recommend our neutralizing buffer.	10 vials of 10ml	
MCNPB10	Neutralizing Peptone Buffer + LTTh		10 vials of 10ml	

MECCONTI

Mecconti S.A.R.L. Sp. z o. o.

Rakowiecka 36 loc.319a, 02-532 Warsaw, Poland

tel. +48 22 397 45 33, fax +48 22 504 72 13 www.mecconti.com, sales@mecconti.com